sebia

CAPILLARYS HEMOGLOBIN(E)

Ref. 2007

CAPILLARYS 2

IVD

CE

2012/06

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PERFORMANCE DATA

NOTE : New CAPILLARYS HEMOGLOBIN(E) performance data have not yet been translated into each language of the instructions for use. They are only available in French and English.

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Precision

The precision of the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP5-A2 guideline "Evaluation of Precision Performance of Quantitative Measurements Methods".

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Seven (7) different blood samples were run using the CAPILLARYS HEMOGLOBIN(E) procedure in all capillaries of the same CAPILLARYS 2 instrument and with 1 lot of CAPILLARYS HEMOGLOBIN(E) kit. The analyzed blood samples included 2 samples with normal Hb A2 level (No. 1 and 5), 2 samples with increased Hb A2 level (No. 2 and 6), 1 sample with low Hb A2 level (No. 3), 1 pathological sample with Hb F and Hb S (No. 4) and 1 sample with increased Hb F level (No. 7). In this study, each blood sample was analyzed on all capillaries from the same instrument, including 40 runs over 20 working days (at 2 different times of the day). Within each run, samples were analyzed in duplicate.

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	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb A %	97.5	96.1	98.2	49.0	97.3	93.5	33.7
Within-run reproducibility (CV %)	0.0	0.0	0.0	0.2	0.0	0.1	0.5
Between-run reproducibility (CV %)	0.0	0.0	0.0	0.1	0.0	0.1	1.1
Between-day reproducibility (CV %)	0.0	0.0	0.0	0.1	0.0	0.1	2.0
Total (CV %)	0.0	0.1	0.1	0.2	0.1	0.1	2.3

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb A2 %	2.5	3.9	1.8	2.8	2.7	6.5	/
Within-run reproducibility (CV %)	1.5	1.0	2.3	1.5	1.7	0.9	/
Between-run reproducibility (CV %)	1.0	0.3	1.5	0.0	1.6	0.7	/
Between-day reproducibility (CV %)	0.0	0.6	0.0	0.4	1.3	0.8	/
Total (CV %)	1.8	1.3	2.7	1.6	2.7	1.5	/

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb F %	/	/	/	12.0	/	/	65.8
Within-run reproducibility (CV %)	/	/	/	0.5	/	/	0.2
Between-run reproducibility (CV %)	/	/	/	0.0	/	/	0.6
Between-day reproducibility (CV %)	/	/	/	0.4	/	/	1.0
Total (CV %)	/	/	/	0.6	/	/	1.1

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb S %	/	/	/	36.2	/	/	/
Within-run reproducibility (CV %)	/	/	/	0.2	/	/	/
Between-run reproducibility (CV %)	/	/	/	0.1	/	/	/
Between-day reproducibility (CV %)	/	/	/	0.2	/	/	/
Total (CV %)	/	/	/	0.3	/	/	/

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The reproducibility study was conducted using 7 different blood samples that were tested on 3 CAPILLARYS 2 instruments with 3 lots of CAPILLARYS HEMOGLOBIN(E) kits. The analyzed blood samples included 2 samples with normal Hb A2 level (No. 1 and 5), 2 samples with increased Hb A2 level (No. 2 and 6), 1 sample with low Hb A2 level (No. 3), 1 pathological sample with Hb F and Hb S (No. 4) and 1 sample with increased Hb F level (No. 7). In this study, each blood sample was analyzed on all capillaries from the 3 CAPILLARYS 2 instruments, including 60 runs over 26 working days (at 2 different times of the day). Within each run, samples were analyzed in duplicate.

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Hemoglobin	Total CV % range
Hb A	0.0 – 3.1
Hb A2	0.0 – 3.7
Hb F	0.3 – 1.5
Hb S	0.3 – 0.6

Linearity

The linearity of the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP6-A guideline "Evaluation of the Linearity of Quantitative Measurement Procedures: A statistical Approach".

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Hb A2 linearity

One Hb A2 enriched blood sample (containing 13.2 g/dL total hemoglobin with 9.9 % Hb A2) was mixed with a Hb A2 depleted blood sample (containing 13.4 g/dL total hemoglobin with 0.0 % Hb A2) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument gave a good linearity for Hb A and Hb A2 within the entire range studied, with a maximum of about 1.3 g/dL for Hb A2 in the analyzed sample (**between 0.0 and 9.9 % of Hb A2**).

Hb F linearity

One umbilical cord blood sample (containing 16.8 g/dL total hemoglobin with 87.8 % Hb F) was mixed with a normal blood sample (containing 10.8 g/dL total hemoglobin with 0.0 % Hb F) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument gave a good linearity for Hb A and Hb F within the entire range studied, with a maximum of about 14.7 g/dL for Hb F in the analyzed sample (**between 0.0 and 87.8 % of Hb F**).

Hb S linearity

One blood sample with Hb S (containing 8.9 g/dL total hemoglobin with 71.3 % Hb S and 0.0 % Hb A) was mixed with a normal blood sample (containing 14.4 g/dL total hemoglobin with 0.0 % Hb S and 97.4 % Hb A) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument gave a good linearity for Hb S within the entire range studied, with a maximum of about 6.3 g/dL for Hb S in the analyzed sample (**between 0.0 and 71.3 % of Hb S**) and a good linearity for Hb A within the entire range studied, with a maximum of about 14.0 g/dL for Hb A in the analyzed sample (**between 0.0 and 97.4 % of Hb A**).

Hb C linearity

One blood sample with Hb C (containing 27.8 g/dL total hemoglobin with 30.0 % Hb C) was mixed with a normal blood sample (containing 25.4 g/dL total hemoglobin with 0.0 % Hb C) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument gave a good linearity for Hb C within the entire range studied, with a maximum of about 8.3 g/dL for Hb C in the analyzed sample (**between 0.0 and 30.0 % of Hb C**).

Hb D linearity

One blood sample with Hb D (containing 24.6 g/dL total hemoglobin with 39.4 % Hb D) was mixed with a normal blood sample (containing 30.5 g/dL total hemoglobin with 0.0 % Hb D) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument gave a good linearity for Hb D within the entire range studied, with a maximum of about 9.7 g/dL for Hb D in the analyzed sample (**between 0.0 and 39.4 % of Hb D**).

Hb E linearity

One blood sample with Hb E (containing 10.8 g/dL total hemoglobin with 95.4 % Hb E) was mixed with a normal blood sample (containing 12.3 g/dL total hemoglobin with 0.0 % Hb E) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 instrument gave a good linearity for Hb E within the entire range studied, with a maximum of about 10.3 g/dL for Hb E in the analyzed sample (**between 0.0 and 95.4 % of Hb E**).

Accuracy – Internal correlation

The internal concordance study of the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP09-A2 guideline "Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition (Interim Revision)".

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The levels of hemoglobin fractions were measured in 56 blood samples, including 20 samples with hemoglobin variants such as hemoglobins S, C, D and E, both by electrophoretic separations obtained with the CAPILLARYS HEMOGLOBIN(E) procedure with the CAPILLARYS 2 instrument and a commercially available HPLC system for hemoglobins analysis.

The measured values of hemoglobin fractions from both procedures were analyzed by a linear regression statistical procedure. The results of linear regression analysis for Hb A, Hb A2, Hb F and Hb S are tabulated below ($y = \text{CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 instrument}$) :

Hb fraction	Number of samples	Correlation coefficient	y-Intercept	Slope	Range of Hb % values CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 instrument
Hb A	56	0.995	- 11.24	1.32	20.9 – 98.1
Hb A2	44	0.977	- 0.07	1.13	0.1 – 6.2
Hb F	56	1.000	- 0.35	0.93	0.0 – 79.0
Hb S	8	0.997	- 0.39	1.07	6.5 – 40.9

All abnormal hemoglobins or abnormal levels of normal hemoglobins detected with the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 instrument were in agreement with the comparative HPLC system. There was no case observed of false positive, i.e., detection of an abnormal band or abnormal level of a normal band where no such abnormality existed.

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In study No. 1, the levels of hemoglobin fractions were measured in 123 blood samples, including 33 samples with hemoglobin variants such as hemoglobins S, C and E, both by electrophoretic separations obtained with the CAPILLARYS HEMOGLOBIN(E) procedure with the CAPILLARYS 2 instrument and a commercially available HPLC system for hemoglobins analysis.

The measured values of hemoglobin fractions from both procedures were analyzed by a linear regression statistical procedure. The results of linear regression analysis for Hb A, Hb A2, Hb F, Hb S and Hb C are tabulated below ($y = \text{CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 instrument}$) :

Hb fraction	Number of samples	Correlation coefficient	y-Intercept	Slope	Range of Hb % values CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 instrument
Hb A	121	0.999	- 1.72	1.14	0.0 – 98.8
Hb A2	93	0.983	- 0.21	1.12	1.0 – 6.6
Hb F	103	0.995	- 0.62	0.99	0.0 – 90.0
Hb S	26	0.988	1.57	0.99	14.7 – 54.9
Hb C	5	0.999	- 0.42	0.99	10.0 – 43.7

All abnormal hemoglobins or abnormal levels of normal hemoglobins detected with the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 instrument were in agreement with the comparative HPLC system. There was no case observed of false positive, i.e., detection of an abnormal band or abnormal level of a normal band where no such abnormality existed.

In study No. 2, the levels of hemoglobin fractions were measured in 183 blood samples, including 83 samples with hemoglobin variants such as hemoglobins S, C and D, both by electrophoretic separations obtained with the CAPILLARYS HEMOGLOBIN(E) procedure with the CAPILLARYS 2 instrument and a commercially available HPLC system for hemoglobins analysis.

The measured values of hemoglobin fractions from both procedures were analyzed by a linear regression statistical procedure. The results of linear regression analysis for Hb A, Hb A2, Hb F, Hb S, Hb C and Hb D are tabulated below ($y = \text{CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 instrument}$) :

Hb fraction	Number of samples	Correlation coefficient	y-Intercept	Slope	Range of Hb % values CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 instrument
Hb A	180	0.991	- 1.75	1.13	0.0 – 98.2
Hb A2	113	0.925	- 0.14	1.02	0.4 – 6.0
Hb F	181	0.991	0.13	1.06	0.0 – 98.0
Hb S	67	0.994	2.49	0.97	7.4 – 93.6
Hb C	13	0.985	1.29	0.92	12.8 – 45.3
Hb D	3	0.984	0.72	1.07	34.5 – 41.1

All abnormal hemoglobins or abnormal levels of normal hemoglobins detected with the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 instrument were in agreement with the comparative HPLC system. There was no case observed of false positive, i.e., detection of an abnormal band or abnormal level of a normal band where no such abnormality existed.

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PERFORMANCE DATA

NOTE : New CAPILLARYS HEMOGLOBIN(E) performance data have not yet been translated into each language of the instructions for use.

They are only available in French and English.

Each translation will be added in the instructions in the next version of the INSTRUCTIONS & SAFETY DATA SHEETS DVD.

Precision

The precision of the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 FLEX-PIERCING instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP5-A2 guideline "Evaluation of Precision Performance of Quantitative Measurements Methods".

The means, standard deviations and coefficients of variation (CV's %) (n = 80) were calculated for percentages (%) of hemoglobin fractions for each sample, using statistical tools recommended by CLSI.

The results obtained with the CAPILLARYS HEMOGLOBIN(E) procedure indicate a very good reproducibility for quantitative analysis for each hemoglobin component. All electrophoregrams were also interpreted visually.

The results presented below have been obtained using the standard parameters of the CAPILLARYS software (smoothing 0 and hemoglobin fractions automatic quantification with HEMOGLOBIN(E) analysis program).

Reproducibility between capillaries from the same instrument

Seven (7) different blood samples were run using the CAPILLARYS HEMOGLOBIN(E) procedure in all capillaries of the same CAPILLARYS 2 FLEX-PIERCING instrument and with 1 lot of CAPILLARYS HEMOGLOBIN(E) kit. The analyzed blood samples included 2 samples with normal Hb A2 level (No. 1 and 5), 2 samples with increased Hb A2 level (No. 2 and 6), 1 sample with low Hb A2 level (No. 3), 1 pathological sample with Hb F and Hb S (No. 4) and 1 sample with increased Hb F level (No. 7). In this study, each blood sample was analyzed on all capillaries from the same instrument, including 40 runs over 20 working days (at 2 different times of the day). Within each run, samples were analyzed in duplicate.

The results for Hb A, Hb A2, Hb F and Hb S percentages are summarized in the following tables.

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb A %	97.2	95.4	98.1	54.9	97.4	93.8	34.3
Within-run reproducibility (CV %)	0.0	0.0	0.0	0.2	0.1	0.1	0.5
Between-run reproducibility (CV %)	0.0	0.0	0.0	0.2	0.0	0.0	0.6
Between-day reproducibility (CV %)	0.0	0.0	0.0	0.1	0.0	0.0	0.9
Total (CV %)	0.0	0.1	0.0	0.4	0.1	0.1	1.2

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb A2 %	2.8	4.7	1.9	2.8	2.6	6.2	/
Within-run reproducibility (CV %)	1.6	1.0	1.5	1.6	2.6	1.3	/
Between-run reproducibility (CV %)	0.4	0.0	1.2	0.4	0.0	0.6	/
Between-day reproducibility (CV %)	0.1	0.7	0.0	0.6	0.7	0.6	/
Total (CV %)	1.6	1.2	1.9	1.7	2.6	1.6	/

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb F %	/	/	/	8.3	/	/	65.2
Within-run reproducibility (CV %)	/	/	/	0.7	/	/	0.3
Between-run reproducibility (CV %)	/	/	/	0.8	/	/	0.4
Between-day reproducibility (CV %)	/	/	/	0.3	/	/	0.5
Total (CV %)	/	/	/	1.1	/	/	0.7

	Sample No. 1	Sample No. 2	Sample No. 3	Sample No. 4	Sample No. 5	Sample No. 6	Sample No. 7
Mean Hb S %	/	/	/	33.9	/	/	/
Within-run reproducibility (CV %)	/	/	/	0.4	/	/	/
Between-run reproducibility (CV %)	/	/	/	0.4	/	/	/
Between-day reproducibility (CV %)	/	/	/	0.1	/	/	/
Total (CV %)	/	/	/	0.6	/	/	/

In addition, none of the repeats showed false positive or false negative values.

Reproducibility between lots and instruments

The reproducibility study was conducted using 7 different blood samples that were tested on 3 CAPILLARYS 2 FLEX-PIERCING instruments with 3 lots of CAPILLARYS HEMOGLOBIN(E) kits. The analyzed blood samples included 2 samples with normal Hb A2 level (No. 1 and 5), 2 samples with increased Hb A2 level (No. 2 and 6), 1 sample with low Hb A2 level (No. 3), 1 pathological sample with Hb F and Hb S (No. 4) and 1 sample with increased Hb F level (No. 7). In this study, each blood sample was analyzed on all capillaries from the 3 CAPILLARYS 2 FLEX-PIERCING instruments, including 60 runs over 27 working days (at 2 different times of the day). Within each run, samples were analyzed in duplicate.

The following table summarizes the total instrument-reagent C.V. % range for the individual hemoglobin Hb A, Hb A2, Hb F and Hb S fractions tested.

Hemoglobin	Total CV % range
Hb A	0.0 – 2.0
Hb A2	0.0 – 5.5
Hb F	0.2 – 1.3
Hb S	0.4 – 0.7

Linearity

The linearity of the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 FLEX-PIERCING instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP6-A guideline "Evaluation of the Linearity of Quantitative Measurement Procedures: A statistical Approach".

The results for percentages (%) of hemoglobin fractions were analyzed using statistical tools recommended by CLSI.

Hb A2 linearity

One Hb A2 enriched blood sample (containing 13.2 g/dL total hemoglobin with 9.9 % Hb A2) was mixed with a Hb A2 depleted blood sample (containing 13.4 g/dL total hemoglobin with 0.0 % Hb A2) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument gave a good linearity for Hb A and Hb A2 within the entire range studied, with a maximum of about 1.3 g/dL for Hb A2 (**between 0.0 and 9.9 % of Hb A2**).

Hb F linearity

One umbilical cord blood sample (containing 12.7 g/dL total hemoglobin with 75.5 % Hb F) was mixed with a normal blood sample (containing 8.8 g/dL total hemoglobin with 0.0 % Hb F) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument gave a good linearity for Hb A and Hb F within the entire range studied, with a maximum of about 9.6 g/dL for Hb F (**between 0.0 and 75.5 % of Hb F**).

Hb S linearity

One blood sample with Hb S (containing 6.0 g/dL total hemoglobin with 90.7 % Hb S and 0.0 % Hb A) was mixed with a normal blood sample (containing 9.1 g/dL total hemoglobin with 0.0 % Hb S and 97.5 % Hb A) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument gave a good linearity for Hb S within the entire range studied, with a maximum of about 5.4 g/dL for Hb S (**between 0.0 and 90.7 % of Hb S**) and a good linearity for Hb A within the entire range studied, with a maximum of about 8.9 g/dL for Hb A (**between 0.0 and 97.5 % of Hb A**).

Hb C linearity

One blood sample with Hb C (containing 11.5 g/dL total hemoglobin with 33.2 % Hb C) was mixed with a normal blood sample (containing 13.8 g/dL total hemoglobin with 0.0 % Hb C) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument gave a good linearity for Hb C within the entire range studied, with a maximum of about 3.8 g/dL for Hb C in the analyzed sample (**between 0.0 and 33.2 % of Hb C**).

Hb D linearity

One blood sample with Hb D (containing 14.3 g/dL total hemoglobin with 40.9 % Hb D) was mixed with a normal blood sample (containing 14.0 g/dL total hemoglobin with 0.0 % Hb D) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument gave a good linearity for Hb D within the entire range studied, with a maximum of about 5.8 g/dL for Hb D in the analyzed sample (**between 0.0 and 40.9 % of Hb D**).

Hb E linearity

One blood sample with Hb E (containing 10.8 g/dL total hemoglobin with 96.1 % Hb E) was mixed with a normal blood sample (containing 12.3 g/dL total hemoglobin with 0.0 % Hb E) within different proportions and the dilutions were electrophoresed with CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument.

The CAPILLARYS HEMOGLOBIN(E) procedure performed with CAPILLARYS 2 FLEX-PIERCING instrument gave a good linearity for Hb E within the entire range studied, with a maximum of about 10.4 g/dL for Hb E in the analyzed sample (**between 0.0 and 96.1 % of Hb E**).

Accuracy – Internal correlation

The internal concordance study of the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 FLEX-PIERCING instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP09-A2 guideline "Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition (Interim Revision)".

The results for percentages (%) of hemoglobin fractions were analyzed using statistical tools recommended by CLSI.

NOTE: The results presented below have been obtained from 1 internal accuracy study that has been performed in SEBIA facility. The analyzed blood samples and their diagnostic assessment were provided by 11 hospital laboratories in France and USA. The diagnosis was based on a routine HPLC procedure.

The levels of hemoglobin fractions were measured in 56 blood samples, including 20 samples with hemoglobin variants such as hemoglobins S, C, D and E, both by electrophoretic separations obtained with the CAPILLARYS HEMOGLOBIN(E) procedure with the CAPILLARYS 2 FLEX-PIERCING instrument and a commercially available HPLC system for hemoglobins analysis.

The measured values of hemoglobin fractions from both procedures were analyzed by a linear regression statistical procedure. The results of linear regression analysis for Hb A, Hb A2, Hb F and Hb S are tabulated below ($y = \text{CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 FLEX-PIERCING instrument}$) :

Hb fraction	Number of samples	Correlation coefficient	y-Intercept	Slope	Range of Hb % values CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 FLEX- PIERCING instrument
Hb A	56	0.996	- 11.42	1.33	21.0 – 98.2
Hb A2	44	0.978	- 0.08	1.13	0.1 – 6.3
Hb F	56	1.000	- 0.38	0.93	0.0 – 79.0
Hb S	8	0.998	0.06	1.07	6.8 – 41.2

All abnormal hemoglobins or abnormal levels of normal hemoglobins detected with the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 FLEX-PIERCING instrument were in agreement with the comparative HPLC system. There was no case observed of false positive, i.e., detection of an abnormal band or abnormal level of a normal band where no such abnormality existed.

Accuracy – External correlations

NOTE: The results presented below have been obtained from 2 external accuracy studies that have been performed in 2 hospital laboratories located in the USA. The diagnosis was based on a routine HPLC procedure.

In study No. 1, the levels of hemoglobin fractions were measured in 123 blood samples, including 33 samples with hemoglobin variants such as hemoglobins S, C and E, both by electrophoretic separations obtained with the CAPILLARYS HEMOGLOBIN(E) procedure with the CAPILLARYS 2 FLEX-PIERCING instrument and a commercially available HPLC system for hemoglobins analysis.

The measured values of hemoglobin fractions from both procedures were analyzed by a linear regression statistical procedure. The results of linear regression analysis for Hb A, Hb A2, Hb F, Hb S and Hb C are tabulated below ($y = \text{CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 FLEX-PIERCING instrument}$) :

Hb fraction	Number of samples	Correlation coefficient	y-Intercept	Slope	Range of Hb % values CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 FLEX- PIERCING instrument
Hb A	121	0.999	- 1.71	1.14	0.0 – 98.8
Hb A2	93	0.986	- 0.31	1.16	0.9 – 6.5
Hb F	103	0.997	- 0.65	0.97	0.0 – 90.3
Hb S	26	0.998	- 0.10	1.05	14.2 – 54.8
Hb C	5	0.999	- 1.08	1.02	9.7 – 44.7

All abnormal hemoglobins or abnormal levels of normal hemoglobins detected with the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 FLEX-PIERCING instrument were in agreement with the comparative HPLC system. There was no case observed of false positive, i.e., detection of an abnormal band or abnormal level of a normal band where no such abnormality existed.

In study No. 2, the levels of hemoglobin fractions were measured in 183 blood samples, including 83 samples with hemoglobin variants such as hemoglobins S, C and D, both by electrophoretic separations obtained with the CAPILLARYS HEMOGLOBIN(E) procedure with the CAPILLARYS 2 FLEX-PIERCING instrument and a commercially available HPLC system for hemoglobins analysis.

The measured values of hemoglobin fractions from both procedures were analyzed by a linear regression statistical procedure. The results of linear regression analysis for Hb A, Hb A2, Hb F, Hb S, Hb C and Hb D are tabulated below ($y = \text{CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 FLEX-PIERCING instrument}$) :

Hb fraction	Number of samples	Correlation coefficient	y-Intercept	Slope	Range of Hb % values CAPILLARYS HEMOGLOBIN(E) with CAPILLARYS 2 FLEX- PIERCING instrument
Hb A	180	0.991	- 1.75	1.13	0.0 – 98.2
Hb A2	113	0.926	- 0.18	1.04	0.4 – 6.0
Hb F	181	0.991	0.11	1.06	0.0 – 98.2
Hb S	67	0.996	0.66	1.00	5.0 – 93.9
Hb C	13	0.996	0.81	0.94	12.5 – 45.7
Hb D	3	0.980	3.04	1.01	34.7 – 40.9

All abnormal hemoglobins or abnormal levels of normal hemoglobins detected with the CAPILLARYS HEMOGLOBIN(E) procedure performed with the CAPILLARYS 2 FLEX-PIERCING instrument were in agreement with the comparative HPLC system. There was no case observed of false positive, i.e., detection of an abnormal band or abnormal level of a normal band where no such abnormality existed.

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Zone	Hémoglobines / Hemoglobins (Hb)
Z1	Hb Santa Ana (chaîne libre alpha), Hb Mizuho (pic mineur), Hb δA'2, Hb αA2, Hb F-Hull, Hb T-Cambodia, variants de Hb A2 "Arya", de Hb A2 "Hasharon", de Hb A2 "Fort de France", de Hb A2 "Ottawa", de Hb A2 "Chad", de Hb A2 "G-Norfolk", de Hb A2 "Matsue-Oki", de Hb A2 "Stanleyville II", de Hb A2 "Montgomery", de Hb A2 "Winnipeg", de Hb A2 "Q-India", de Hb A2 « G-Pest", de Hb A2 "Memphis", de Hb A2 "Inkster", de Hb A2 "Chapel Hill" et de Hb A2 "Q-Thailand"
	<i>Hb Santa Ana (free alpha chain), Hb Mizuho (minor peak), Hb δA'2, Hb αA2, Hb F-Hull, Hb T-Cambodia, "Arya" Hb A2 variant, "Hasharon" Hb A2 variant, "Fort de France" Hb A2 variant, "Ottawa" Hb A2 variant, "Chad" Hb A2 variant, "G-Norfolk" Hb A2 variant, "Matsue-Oki" Hb A2 variant, "Stanleyville II" Hb A2 variant, "Montgomery" Hb A2 variant, "Winnipeg" Hb A2 variant, "Q-India" Hb A2 variant, "G-Pest" Hb A2 variant, "Memphis" Hb A2 variant, "Inkster" Hb A2 variant, "Chapel Hill" Hb A2 variant, "Q-Thailand" Hb A2 variant</i>
Z2	Hb C, Hb F-Texas, Hb Constant Spring, Hb C-Harlem (C-Georgetown), variant de Hb A2 "Setif", variant de Hb A2 "Bassett", variant de Hb A2 "Swan River", variant de Hb A2 "Manitoba I", variant de Hb A2 "Manitoba II"
	<i>Hb C, Hb F-Texas, Hb Constant Spring, Hb C-Harlem (C-Georgetown), "Setif" Hb A2 variant, "Bassett" Hb A2 variant, "Swan River" Hb A2 variant, "Manitoba I" Hb A2 variant, "Manitoba II" Hb A2 variant</i>
Z3	Hb A2, Hb O-Arab, Hb Chad (E-Keelung), Hb E-Saskatoon
	<i>Hb A2, Hb O-Arab, Hb Chad (E-Keelung), Hb E-Saskatoon</i>
Z4	Hb E, Hb Köln (Ube-1), Hb Buenos Aires (pic mineur), Hb Agenogi, Hb G-Siriraj, Hb Santa Ana, Hb A2-Babinga, Hb M-Saskatoon (pic mineur), variant de Hb A2 "M-Iwate", Hb C dégradée
	<i>Hb E, Hb Köln (Ube-1), Hb Buenos Aires (minor peak), Hb Agenogi, Hb G-Siriraj, Hb Santa Ana, Hb A2-Babinga, Hb M-Saskatoon (minor peak), "M-Iwate" Hb A2 variant, denatured Hb C</i>
Z5	Hb S, Hb Dhofar (Yukuhashi), Hb Arya, Hb Hasharon (Sinai), Hb Handsworth, Hb Ottawa (Siam), Hb S-Antilles, Hb Fort de France, Hb Hamadan, Hb Montgomery, variant de Hb A2 "Lombard", variant de Hb A2 "Cemenelum", variant de Hb A2 "Jackson", Hb O-Arab dégradée
	<i>Hb S, Hb Dhofar (Yukuhashi), Hb Arya, Hb Hasharon (Sinai), Hb Handsworth, Hb Ottawa (Siam), Hb S-Antilles, Hb Fort de France, Hb Hamadan, Hb Montgomery, "Lombard" Hb A2 variant, "Cemenelum" Hb A2 variant, "Jackson" Hb A2 variant, denatured Hb O-Arab</i>
Z6	Hb D-Punjab (D-Los Angeles), Hb Stanleyville II, Hb Osu Christiansborg, Hb Leiden, Hb G-Philadelphie, Hb Muravera, Hb D-Bushman, Hb G-Norfolk, Hb Matsue-Oki, Hb D-Ouled Rabah, Hb Lepore (Lepore-BW), Hb Muskegen, Hb D-Ibadan, Hb Buenos Aires (pic mineur), Hb Summer Hill, Hb Q-India, Hb Q-Iran, Hb Korle-Bu (G-Accra), Hb Köln, Hb Fort Worth, Hb G-Taïpei, Hb Winnipeg, Hb G-Siriraj, Hb G-Pest, Hb D-Iran, Hb G-Coushatta (G-Saskatoon), Hb Inkster, Hb Setif, Hb Memphis, Hb P-Nilotic, Hb Willamette, variant de Hb A2 "J-Rajappen", variant de Hb A2 "J-Anatolia", variant de Hb A2 "J-Oxford", variant de Hb A2 "J-Broussais", variant de Hb A2 "J-Toronto", variant de Hb A2 "Mexico", variant de Hb A2 "J-Habana", variant de Hb A2 "J-Rovigo", Hb E dégradée
	<i>Hb D-Punjab (D-Los Angeles), Hb Stanleyville II, Hb Osu Christiansborg, Hb Leiden, Hb G-Philadelphia, Hb Muravera, Hb D-Bushman, Hb G-Norfolk, Hb Matsue-Oki, Hb D-Ouled Rabah, Hb Lepore (Lepore-BW), Hb Muskegen, Hb D-Ibadan, Hb Buenos Aires (minor peak), Hb Summer Hill, Hb Q-India, Hb Q-Iran, Hb Korle-Bu (G-Accra), Hb Köln, Hb Fort Worth, Hb G-Taïpei, Hb Winnipeg, Hb G-Siriraj, Hb G-Pest, Hb D-Iran, Hb G-Coushatta (G-Saskatoon), Hb Inkster, Hb Setif, Hb Memphis, Hb P-Nilotic, Hb Willamette, "J-Rajappen" Hb A2 variant, "J-Anatolia" Hb A2 variant, "J-Oxford" Hb A2 variant, "J-Broussais" Hb A2 variant, "J-Toronto" Hb A2 variant, "Mexico" Hb A2 variant, "J-Habana" Hb A2 variant, "J-Rovigo" Hb A2 variant, denatured Hb E</i>
Z7	Hb F, Hb Q-Thailand (G-Taichung), Hb Alabama, Hb Chapel Hill, Hb Bassett, Hb G-San José, Hb Richmond, Hb Barcelona, Hb Geldrop Santa Anna, Hb Porto Alegre, Hb Swan River, Hb Presbyterian, Hb Burke, Hb Manitoba I, Hb Manitoba II, variant de Hb A2 "J-Paris-I", Hb S dégradée
	<i>Hb F, Hb Q-Thailand (G-Taichung), Hb Alabama, Hb Chapel Hill, Hb Bassett, Hb G-San Jose, Hb Richmond, Hb Barcelona, Hb Geldrop Santa Anna, Hb Porto Alegre, Hb Swan River, Hb Presbyterian, Hb Burke, Hb Manitoba I, Hb Manitoba II, "J-Paris-I" Hb A2 variant, denatured Hb S</i>
Z8	Hb F acétylée, Hb Hinsdale, Hb Alberta, Hb Kempsey, Hb Atlanta, Hb Athens-GA (Waco)
	<i>Acetylated Hb F, Hb Hinsdale, Hb Alberta, Hb Kempsey, Hb Atlanta, Hb Athens-GA (Waco)</i>
Z9	Hb A, Hb Gorwihl (Hinchinbrooke), Hb Phnom Penh, Hb Silver Springs, Hb La Coruna, Hb Bougardirey-Mali, Hb Austin, Hb Buenos Aires (pic majeur), Hb Chicago, Hb Toulon, Hb Okayama, Hb Fontainebleau, Hb Raleigh, Hb Hekinan, Hb Mosella, Hb Dallas, Hb Aztec, Hb Little Rock, Hb Frankfurt, Hb Bethesda, Hb M-Boston (M-Osaka), Hb Brisbane (Great Lakes), Hb Mizuho, Hb Grange Blanche, Hb San Diego, Hb M-Saskatoon (pic majeur), Hb Malmö, Hb Minneapolis Laos, Hb Syracuse, Hb Camperdown
	<i>Hb A, Hb Gorwihl (Hinchinbrooke), Hb Phnom Penh, Hb Silver Springs, Hb La Coruna, Hb Bougardirey-Mali, Hb Austin, Hb Buenos Aires (major peak), Hb Chicago, Hb Toulon, Hb Okayama, Hb Fontainebleau, Hb Raleigh, Hb Hekinan, Hb Mosella, Hb Dallas, Hb Aztec, Hb Little Rock, Hb Frankfurt, Hb Bethesda, Hb M-Boston (M-Osaka), Hb Brisbane (Great Lakes), Hb Mizuho, Hb Grange Blanche, Hb San Diego, Hb M-Saskatoon (major peak), Hb Malmö, Hb Minneapolis Laos, Hb Syracuse, Hb Camperdown</i>
Z10	Hb Hope, Hb M-Iwate (M-Kankakee), Hb Camden (Tokuchi)
	<i>Hb Hope, Hb M-Iwate (M-Kankakee), Hb Camden (Tokuchi)</i>

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
Zone	Hémoglobines / Hemoglobins (Hb)
Z11	Hb A dégradée, Hb Providence (pic X-Asn), Hb K-Woolwich, Hb Lombard, Hb Kaohsiung (New York), Hb Fannin Lubbock, Hb Andrew Minneapolis, Hb Jackson, Hb Himeji, variant de Hb A2 "I (I-Texas)"
	<i>Denatured Hb A, Hb Providence (X-Asn peak), Hb K-Woolwich, Hb Lombard, Hb Kaohsiung (New York), Hb Fannin Lubbock, Hb Andrew Minneapolis, Hb Jackson, Hb Himeji, "I (I-Texas)" Hb A2 variant</i>
Z12	Hb Bart, Hb Cemenelum, Hb J-Calabria (J-Bari), Hb Providence (pic X-Asp), Hb J-Rajappen, Hb Grady, Hb J-Anatolia, Hb J-Broussais (Tagawa-I), Hb J-Chicago, Hb J-Oxford (I-Interlaken), Hb J-Toronto, Hb J-Meinung (J-Bangkok), Hb Ube-2, Hb Mexico (J-Paris-II), Hb J-Habana, Hb J-Baltimore (N-New Haven), Hb J-Paris-I (J-Aljezur)
	<i>Hb Bart's, Hb Cemenelum, Hb J-Calabria (J-Bari), Hb Providence (X-Asp peak), Hb J-Rajappen, Hb Grady, Hb J-Anatolia, Hb J-Broussais (Tagawa-I), Hb J-Chicago, Hb J-Oxford (I-Interlaken), Hb J-Toronto, Hb J-Meinung (J-Bangkok), Hb Ube-2, Hb Mexico (J-Paris-II), Hb J-Habana, Hb J-Baltimore (N-New Haven), Hb J-Paris-I (J-Aljezur)</i>
Z13	Hb J-Rovigo, Hb N-Baltimore (Hopkins-I), Hb J-Norfolk (Kagoshima), Hb J-Kaohsiung (J-Honolulu)
	<i>Hb J-Rovigo, Hb N-Baltimore (Hopkins-I), Hb J-Norfolk (Kagoshima), Hb J-Kaohsiung (J-Honolulu)</i>
Z14	Hb N-Seattle
	<i>Hb N-Seattle</i>
Z15	Hb H, Hb I (I-Texas)
	<i>Hb H, Hb I (I-Texas)</i>

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








Figure 9


Sang avec variant Hb D-Punjab
Blood sample with Hb D-Punjab variant

Figure 10

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







Zone	Hémoglobines / Hemoglobins (Hb)
C1	Hb C
	<i>Hb C</i>
C2	Hb A2, Hb O-Arab
	<i>Hb A2, Hb O-Arab</i>
C3	Hb E, méthémoglobine, Hb C dégradée
	<i>Hb E, methemoglobin, denatured Hb C</i>
C4	Hb S
	<i>Hb S</i>
C5	Hb D-Punjab (D-Los Angeles), Hb G-Philadelphie, Hb Korle-Bu (G-Accra), Hb E dégradée
	<i>Hb D-Punjab (D-Los Angeles), Hb G-Philadelphia, Hb Korle-Bu (G-Accra), denatured Hb E</i>
C6	Hb S dégradée
	<i>Denatured Hb S</i>
C7	Hb F
	<i>Hb F</i>
C9	Hb F dégradée
	<i>Denatured Hb F</i>
C10	Hb A
	<i>Hb A</i>
C11	Hb A dégradée, Hb Hope
	<i>Denatured Hb A, Hb Hope</i>
C12	Hb Bart
	<i>Hb Bart's</i>

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